

KFT-41V Karl Fischer Titrator

GENERAL FEATURES

- Color High contrast LCD touchscreen, 7.0 inches.
- Coulometric KF Titration.
- Support oven sampling method, to meet the gas, solid, liquid in the constant moisture, trace moisture measurement.

Volumetric Karl Fischer Titration

- Support Auto Titration Method, KF Titer Determination Method, Oven Titration Method and Oven Blank Titration Method.
- The solvent manager adopts anti-leakage liquid and anti-pour design of waste liquid bottle, supports to fill titrant, drain solution, feed solvent, mix titrant and clean burette, ensure safe handling of Karl Fischer chemicals.

Includes

- Solvent manager
- 10mL burette
- Volumetric vessel
- Measuring electrode



Data management

- Auto/manual drift adjust ensures accurate results.
- Selectable units including μg , mg, %, ppm, mg/L and $\mu\text{g/mL}$.
- User password management.
- Store up to 200 titration data sets (GLP-compliant).
- Support for RS-232 communication for print and the print format can be alternate.
- Support for USB communication for PC software.
- Support for U disk for data export with CSV format.
- It needs GMP software to implement GMP management functions.

Model		KFT-41V
Principle		Karl Fischer Volumetric Titration
Water	Range	(0.1~250.0)mg
	Resolution	1 μg
	Unit	μg , mg, %, ppm, $\mu\text{g/mL}$, mg/L
	Repeatability	$\leq 0.5\%$
mV Range		(0~2000)mV
mV Resolution		1mV
Polarization Current Accuracy		$\pm 3\%$
Polarization Current Fluctuation		$\pm 2.5\%/30\text{min}$
Working Current Accuracy		/
Working Current Fluctuation		/
Power Supply		(100~240)VAC, (47~63)Hz
Dimensions		240×370×270mm
Net Weight		About 4kg